

Product Datasheet

p65 (phospho-T60) peptide (orb304903)

Description	The peptide is used to p65 (phospho-T60) antibody ractivity
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	WNK1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P83741 , Q9JIH7 , Q9H4A3
Entrez	232341 , 65125 , 116477
Expiration Date	12 months from date of receipt.