

Product Datasheet

Alpha-adducin (phospho-T445) peptide (orb304918)

Description	The peptide is used to Alpha-adducin (phospho-T445) antibody ractivity
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	ADD1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q63028 , P35611 , Q9QYC0
Entrez	11518 , 118 , 24170
Expiration Date	12 months from date of receipt.