

Product Datasheet

AKT (phospho-Y474) peptide (orb304921)

Description	The peptide is used to AKT (phospho-Y474) antibody ractivity
Reactivity	Human, Mouse, Porcine, Primate, Rat, Sheep, Zebrafish
Tested Applications	Blocking
Target	AKT1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P47196 , P31750 , P31749
Entrez	207 , 24185 , 11651
Expiration Date	12 months from date of receipt.