

## Product Datasheet

### Ephrin B1/2 peptide (orb304983)

<b>Description</b>	The peptide is used to Ephrin B1/2 antibody ractivity
<b>Reactivity</b>	Human, Mouse, Rat, Zebrafish
<b>Tested Applications</b>	Blocking
<b>Target</b>	EFNB1; EFNB2
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Shipped at 4°C. Store at -20°C for one year.
<b>Note</b>	For research use only
<b>Application notes</b>	Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
<b>Purity</b>	> 85%
<b>Source</b>	Synthetic
<b>Uniprot ID</b>	<a href="#">P52796</a> , <a href="#">P52800</a> , <a href="#">P52799</a> , <a href="#">P98172</a> , <a href="#">P52795</a>
<b>Entrez</b>	<a href="#">1948</a> , <a href="#">13641</a> , <a href="#">13642</a> , <a href="#">1947</a>
<b>Expiration Date</b>	12 months from date of receipt.