

## Product Datasheet

### 4EBP1 (phospho-T69) peptide (orb304989)

**Description** The peptide is used to 4EBP1 (phospho-T69) antibody ractivity

**Reactivity** Human, Mouse, Primate, Rat

**Tested Applications** Blocking

**Target** EIF4EBP1

**Form/Appearance** Lyophilized powder

**Storage** Shipped at 4°C. Store at -20°C for one year.

**Note** For research use only

**Application notes** Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.

**Purity** > 85%

**Source** Synthetic

**Uniprot ID** [Q13541](#), [Q62622](#), [Q60876](#)

**Entrez** [13685](#), [116636](#), [1978](#)

**Expiration Date** 12 months from date of receipt.