



Product Datasheet

PIASy1 antibody (orb30521)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	PIAS4
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Application notes	The data on the recent batch and one customer feedback were not positive for Western Blot. We would be able to offer trial sizes at stock at a discount to customers for further testing.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	56504
Uniprot ID	Q8N2W9
NCBI	NP_056981.2
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of bacterial cell
...