




Product Datasheet

PPT1 antibody (orb305764)

Description	Mouse monoclonal antibody to PPT1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Recombinant Protein
Target	PPT1
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	34193
Uniprot ID	P50897
Dilution Range	WB: 1:1000, IHC-P: 1:25, FACS: 1:25
Expiration Date	12 months from date of receipt.



IHC-P staining of human cerebellum secti...



Immunohistochemical staining of mouse ce...