




Product Datasheet

RBX1 antibody (orb305784)

Description	Rabbit polyclonal antibody to RBX1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	RBX1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	12274
Uniprot ID	P62877
Dilution Range	IF/ICC: 1:25, IHC-P: 1:25, WB: 1:1000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse te...



Immunohistochemical staining of paraffin...