




Product Datasheet

AK4 antibody (orb305785)

Description	Rabbit polyclonal antibody to AK4
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	AK4
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	25268
Uniprot ID	P27144
Dilution Range	IF/ICC: 1:25, IHC-P: 1:25, WB: 1:1000
Expiration Date	12 months from date of receipt.



IHC-P staining of mouse kidney section t...



Immunohistochemical staining of human ki...