




Product Datasheet

Lysozyme antibody (orb305834)

Description	Rabbit polyclonal antibody to Lysozyme
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Synthetic Peptide
Target	LYZ
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	16537
Uniprot ID	P61626
Dilution Range	IHC-P: 1:25, WB: 1:1000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat smal...



Immunohistochemical staining of paraffin...