




Product Datasheet

RAC1 antibody (orb305856)

Description	Mouse monoclonal antibody to RAC1
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Target	RAC1
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	21450
Uniprot ID	P63000
Dilution Range	IHC-P: 1:25, FACS: 1:100, WB: 1:1000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human sk...



Immunohistochemical staining of paraffin...