

Product Datasheet

FABP3 antibody (orb305858)

Description Rabbit polyclonal antibody to FABP3

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications FC, IHC-P, WB

Target FABP3

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles

Note For research use only

Isotype Rabbit IgG



Clonality Polyclonal

MW 14858


Uniprot ID [P05413](#)

Dilution Range IHC-P: 1:25, FACS: 1:25, WB: 1:1000

Expiration Date 12 months from date of receipt.



Immunohistochemical staining of human he...



Immunohistochemical staining of paraffin...