

Product Datasheet

Ribosome-inactivating protein gynostemmin antibody (FITC) (orb30598)

Description	Ribosome-inactivating protein gynostemmin antibody (FITC)
Species/Host	Rabbit
Reactivity	Plant
Conjugation	FITC
Immunogen	Recombinant Gynostemma pentaphyllum Ribosome-inactivating protein gynostemmin protein (24-300AA)
Target	Ribosome-inactivating protein gynostemmin
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5F0I3
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.