

Product Datasheet

MRGPRX3 antibody (orb306076)

Description	Rabbit polyclonal antibody to MRGPRX3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	MRGPRX3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	36483
Uniprot ID	Q96LB0
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

