



Product Datasheet

SULT1A1 antibody (orb306308)

Description	Rabbit polyclonal antibody to SULT1A1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Target	SULT1A1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	34165
Uniprot ID	P50225
Dilution Range	FACS: 1:25, WB: 1:1000
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of HepG2 cells

...