

Product Datasheet

sh3bp4a antibody (orb306338)

Description	Rabbit polyclonal antibody to sh3bp4a
Species/Host	Rabbit
Reactivity	Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Target	sh3bp4a
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	108453
Uniprot ID	Q1LVQ2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

