



Product Datasheet

CAPN6 antibody (orb306343)

Description	Mouse monoclonal antibody to CAPN6
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Target	CAPN6
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	74576
Uniprot ID	Q9Y6Q1
Dilution Range	FACS: 1:25, WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of human placenta ...