




## Product Datasheet

### Srms antibody (orb306379)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Rabbit polyclonal antibody to Srms   |
| <b>Species/Host</b>        | Rabbit   |
| <b>Reactivity</b>          | Mouse  |
| <b>Conjugation</b>         | Unconjugated   |
| <b>Tested Applications</b> | IF, WB   |
| <b>Target</b>              | Srms   |
| <b>Form/Appearance</b>     | Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification. |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles                                       |
| <b>Note</b>                | For research use only  |
| <b>Isotype</b>             | Rabbit IgG   |
| <b>Clonality</b>           | Polyclonal   |
| <b>MW</b>                  | 55758  |
| <b>Uniprot ID</b>          | <a href="#">Q62270</a>   |
| <b>Dilution Range</b>      | WB: 1:2000, WB: 1:1000   |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |



Western blot analysis of mouse liver, mo...



Immunofluorescence analysis of mouse ski...