



Product Datasheet

FGF4 antibody (orb306403)

Description	Rabbit polyclonal antibody to FGF4
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	FGF4
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	22048
Uniprot ID	P08620
Dilution Range	WB: 1:4000
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cell line

...