




Product Datasheet

ATG4A antibody (orb306410)

Description	Mouse monoclonal antibody to ATG4A
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Target	ATG4A
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	IgG2b,k
Clonality	Monoclonal
MW	45378
Uniprot ID	Q8WYN0
Dilution Range	IF/ICC: 1:25, IHC-P: 1:25, FACS: 1:25, WB: 1:500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human br...



Flow cytometric analysis of HeLa cells u...