




Product Datasheet

FES antibody (orb306445)

Description	Mouse monoclonal antibody to FES
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Target	FES
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	IgG1,k
Clonality	Monoclonal
MW	93497
Uniprot ID	P07332
Dilution Range	WB: 1:2000, WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of U87-MG whole ce...



Immunohistochemical staining of human sk...