




Product Datasheet

FAM168B antibody (orb306456)

Description	Rabbit polyclonal antibody to FAM168B
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Target	FAM168B
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	20324
Uniprot ID	A1KXE4
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human br...



Flow cytometric analysis of SH-SY5Y cell...