




## Product Datasheet

### ABCD2 antibody (orb30646)

|                            |                                                                                                                                                                              |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Peptide affinity purified rabbit polyclonal antibody.                                                                                                                        |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                       |
| <b>Reactivity</b>          | Human                                                                                                                                                                        |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                 |
| <b>Tested Applications</b> | FC, IHC-P, WB                                                                                                                                                                |
| <b>Immunogen</b>           | Synthetic Peptide                                                                                                                                                            |
| <b>Target</b>              | ABCD2                                                                                                                                                                        |
| <b>Form/Appearance</b>     | Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification. |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles                                       |
| <b>Note</b>                | For research use only                                                                                                                                                        |
| <b>Isotype</b>             | Rabbit IgG                                                                                                                                                                   |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                   |
| <b>MW</b>                  | 83233                                                                                                                                                                        |
| <b>Uniprot ID</b>          | <a href="#">Q9UBJ2</a>                                                                                                                                                       |
| <b>NCBI</b>                | <a href="#">NP_005155.1</a>                                                                                                                                                  |
| <b>Dilution Range</b>      | WB: 1:1000, IHC-P: 1:10-50, FACS: 1:10-50                                                                                                                                    |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                              |



Flow cytometric analysis of ABCD2 antibody.



Immunohistochemical analysis of formalin-fixed tissue.