




Product Datasheet

CSF1R antibody (orb306507)

Description	Mouse monoclonal antibody to CSF1R
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Target	CSF1R
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	IgG1,k
Clonality	Monoclonal
MW	107984
Uniprot ID	P07333
Dilution Range	IHC-P: 1:25, FACS: 1:25, WB: 1:4000
Expiration Date	12 months from date of receipt.



Western blot analysis of U-87MG whole ce...



Immunohistochemical staining of human sk...