




Product Datasheet

Med12 antibody (orb306655)

Description	Rabbit polyclonal antibody to Med12
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Target	Med12
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	244561
Uniprot ID	A2AGH6
Dilution Range	WB: 1:2000, WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of NIH/3T3 whole c...



Western blot analysis of All lanes using ...