




Product Datasheet

PPP2R1B antibody (orb306717)

Description	Mouse monoclonal antibody to PPP2R1B
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Target	PPP2R1B
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	IgG1,K
Clonality	Monoclonal
MW	66214
Uniprot ID	P30154
Dilution Range	WB: 1:4000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human sp...



Flow cytometric analysis of Jurkat cells...