


Product Datasheet

HRPT2/CDC73 Antibody (orb308827)

Description	HRPT2/CDC73 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	E.coli-derived human HRPT2 recombinant protein (Position: Q119-R520). Human HRPT2 shares 100% and 99.8% amino acid (aa) sequence identity with mouse and rat HRPT2, respectively.
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	60577 MW
Uniprot ID	Q6P1J9
Expiration Date	12 months from date of receipt.



WB analysis of HRPT2 using anti-HRPT2 an...