



Product Datasheet

DUT antibody (orb308875)


Description	Rabbit polyclonal antibody to DUT
Species/Host	Goat
Reactivity	Bovine, Human, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	dUTPase / DUT
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 20kDa band observed in lysates of cell lines HEK293 and HeLa (calculated MW of 17.7kDa according to NP_001939.1). Recommended concentration: 0.1-0.3µg/ml. Primary incubation was 1 hour.. IHC: Paraffin embedded Human Testis and Tonsil. Recommended concentration: 5µg/ml.
Protein Sequence	QLRFARLSEHATAPT
Clonality	Polyclonal
MW	26.6; 17.7; 15.4
Entrez	1854
Dilution Range	ELISA: 1:16000, WB: 0.1-0.3 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293 (Lane 1)...



Immunohistochemical staining of Human Testis...



Immunohistochemical staining of Human Tonsil...