

Product Datasheet

Bcl-6 (Phospho-Ser333) Polyclonal Antibody (orb308962)

| | |
|----------------------------|-------------------------------------------------------------------------------------------|
| Description | Bcl-6 (Phospho-Ser333) Polyclonal Antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, WB |
| Storage | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |
| Note | For research use only |
| Clonality | Recombinant |
| Uniprot ID | P41182 |
| Expiration Date | 12 months from date of receipt. |