

Product Datasheet

SH-PTP2 (Phospho-Tyr542) Polyclonal Antibody (orb308981)

Description	SH-PTP2 (Phospho-Tyr542) Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q06124
Expiration Date	12 months from date of receipt.