

Product Datasheet

Lyn (Phospho-Tyr508) Polyclonal Antibody (orb308990)

Description	Lyn (Phospho-Tyr508) Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P07948
Expiration Date	12 months from date of receipt.