

Product Datasheet

PLC β 3 (Phospho-Ser1105) Polyclonal Antibody (orb308991)

Description	PLC β 3 (Phospho-Ser1105) Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q01970
Expiration Date	12 months from date of receipt.