

Product Datasheet

Histone H1 (Phospho-Thr17) Polyclonal Antibody (orb309066)

Description	Histone H1 (Phospho-Thr17) Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Monkey, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P16401
Expiration Date	12 months from date of receipt.