

Product Datasheet

Histone H2A (Phospho-Thr121) Polyclonal Antibody (orb309082)

| | |
|----------------------------|---|
| Description | Histone H2A (Phospho-Thr121) Polyclonal Antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IHC, WB |
| Storage | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |
| Note | For research use only |
| Clonality | Recombinant |
| Uniprot ID | POC0S8 |
| Expiration Date | 12 months from date of receipt. |