

Product Datasheet

Cleaved-Caspase-4 p20 (Q81) Polyclonal Antibody (orb310049)

Description	Cleaved-Caspase-4 p20 (Q81) Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P49662
Expiration Date	12 months from date of receipt.