

Product Datasheet

GST Tag antibody (orb31021)

Description	Rabbit polyclonal antibody to Bax is a member of the Bcl-2 family and is expressed in a wide variet...
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	This GST Tag antibody is generated from rabbits immunized with a recombinant full length of GST.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

