



Product Datasheet

Kv3.4 Polyclonal Antibody (orb310639)

Description	Kv3.4 Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Monkey, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q03721
Expiration Date	12 months from date of receipt.



Western Blot analysis of SH-SY5Y cells u...



Western Blot analysis of Jurkat cells us...