



Product Datasheet

Bmf BH3 Domain antibody (orb31143)

Description	Rabbit polyclonal antibody to RAD9 is a cell-cycle checkpoint protein with 3-5 exonuclease activity...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	BMF
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	20508
Uniprot ID	Q96LC9
NCBI	NP_277038.1 , NP_001003940.1 , NP_001003942.1 , NP_001003943.1
Dilution Range	WB: 1:1000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...