

## Product Datasheet

### Mouse Acrp30 protein (C-His) (orb311646)

<b>Description</b>	Recombinant Mouse Adipocyte complement-related 30 kDa protein is produced by our Mammalian expressio...
<b>Endotoxins</b>	Less than 0.1 ng/ug (1 EU/ug)
<b>Form/Appearance</b>	Lyophilized from a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
<b>Storage</b>	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
<b>Tag</b>	C-His
<b>Note</b>	For research use only
<b>Application notes</b>	6His tag at the C-terminus.
<b>Uniprot ID</b>	<a href="#">Q60994</a>
<b>Expiration Date</b>	6 months from date of receipt.