

Product Datasheet

Human ABP1 protein (C-His) (orb311651)

Description	Recombinant Human Amiloride-binding protein 1 is produced by our Mammalian expression system and the...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM TrisHCl, 150mM NaCl, 10% Glycerol, pH 7.5.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Uniprot ID	AAH14093.1
Expiration Date	6 months from date of receipt.