

Product Datasheet

Human DUSP3 protein (N-His) (orb311734)

Description	Recombinant Human Vaccinia Virus VH1-related Phosphatase is produced by our E.coli expression system...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Supplied as a 0.2 um filtered solution of PBS, pH 7.4.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-His
Note	For research use only
Application notes	6His tag at the N-terminus.
Uniprot ID	P51452
Expiration Date	6 months from date of receipt.