

Product Datasheet

Human Cell Proliferation-Inducing Gene 15 protein (N-His) (orb311748)

Description	Recombinant Human Ferritin heavy chain is produced by our E.coli expression system and the target ge...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM Tris HCl,150mM NaCl,1mM EDTA,pH8.0 .
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	N-His
Note	For research use only
Application notes	6His tag at the N-terminus.
Uniprot ID	P02794
Expiration Date	6 months from date of receipt.