

Product Datasheet

Human LIPF protein (C-His) (orb311771)

Description	Recombinant Human Gastric Triacylglycerol Lipase is produced by our Mammalian expression system and ...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Supplied as a 0.2 um filtered solution of 25mM Tris-HCl, 100mM glycine, 10% Glycerol, pH 7.3.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Uniprot ID	AAI12273.1
Expiration Date	6 months from date of receipt.