

Product Datasheet

Human GDNF Family Receptor Alpha-2 protein (C-Fc-His) (orb311772)

Description	Recombinant Human Glial Cell line-derived Neurotrophic Factor Receptor alpha 2 is produced by our Ma...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 8.0.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-Fc&His
Note	For research use only
Application notes	Fc, 6His tag at the C-terminus.
Uniprot ID	O00451
Expiration Date	6 months from date of receipt.