

Product Datasheet

Human High Mobility Group Protein B1 protein (orb311803)

Description	Recombinant Human High mobility group protein B1 is produced by our E.coli expression system and the...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 50mM HEPES-Na, 500mM NaCl, 0.6mM DTT, p H 7.9.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	N-MARI
Note	For research use only
Application notes	MARI tag at the N-terminus.
Uniprot ID	P09429
Expiration Date	6 months from date of receipt.