

Product Datasheet

Human IL-1R3 protein (C-Fc-His) (orb311810)

Description	Recombinant Human Interleukin-1 Receptor Accessory Protein is produced by our Mammalian expression s...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-Fc&His
Note	For research use only
Application notes	Fc, 6His tag at the C-terminus.
Uniprot ID	Q9NPH3
Expiration Date	6 months from date of receipt.