

Product Datasheet

Human IDH1 protein (C-His) (orb311859)

Description	Recombinant Human Isocitrate Dehydrogenase [NADP] Cytoplasmic is produced by our Mammalian expressio...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Supplied as a 0.2 um filtered solution of PBS, pH 7.4.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Uniprot ID	O75874
Expiration Date	6 months from date of receipt.