

Product Datasheet

Human Kelch-Like Protein 41 protein (C-His) (orb311861)

| | |
|--------------------------|--|
| Description | Recombinant Human Kelch-like protein 41 is produced by our Mammalian expression system and the target... |
| Endotoxins | Less than 0.1 ng/ug (1 EU/ug) |
| Form/Appearance | Supplied as a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 7.4. |
| Storage | Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles. |
| Tag | C-His |
| Note | For research use only |
| Application notes | 6His tag at the C-terminus. |
| Uniprot ID | P23276 |
| Expiration Date | 6 months from date of receipt. |