

Product Datasheet

Human LTA4H protein (C-His) (orb311872)

Description	Recombinant Human Leukotriene A-4 hydrolase is produced by our Mammalian expression system and the t...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM MES, 2mM EDTA, 20% Glycerol, pH 5.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Uniprot ID	P09960
Expiration Date	6 months from date of receipt.