

Product Datasheet

Human Peptidyl-prolyl Cis-Trans Isomerase FKBP3 protein (N-His) (orb311920)

Description	Recombinant Human 25 kDa FK506-binding protein is produced by our E.coli expression system and the t...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 1mM DTT, 10% Glycerol, pH 8.0 .
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-His
Note	For research use only
Application notes	6His tag at the N-terminus.
Uniprot ID	Q00688
Expiration Date	6 months from date of receipt.